



## **BIOLOGICAL SAFETY CABINETS**


**Maureen O'Leary, PhD, MBA, CBSP**  
**Assistant Director Academic Safety &**  
**Environmental Health**

## BIOSAFETY CABINETS


- Developed for working safely with infectious materials
- HEPA filter=High Efficiency Particulate Air 99.97% min. particle removal for .3microns
- Laminar flow hoods (Clean Benches) -for plant tissue culture, media preparation ONLY
- Fume hoods -for work with volatile chemical compounds




# BIOSAFETY CABINET




## LAMINAR FLOW OR CLEAN BENCH HORIZONTAL FLOW





# LAMINAR FLOW OR CLEAN BENCH VERTICAL FLOW



# BIOSAFETY CABINET AIRFLOW DIAGRAM



# LAMINAR FLOW OR CLEAN BENCHES AIRFLOW DIAGRAM



## BIOSAFETY CABINET TYPES

- Class I: No product protection, equipment use only
- Class II Type A: Both personnel & product protection
  - not suitable for use of volatile and toxic chemicals
  - 70% recirculated air
- Class II Type B1: Microbes plus low level toxic chemicals
  - 70% exhausted, 30% recirculated
- Class II Type B2: no recirculated air
- Class III: Totally enclosed





# BIOSAFETY CABINETS

## PROCEDURES FOR USE

- Turn on cabinet for 5 minutes before initiation
- Disinfect surfaces
- Assemble & organize material clean-contaminated areas, equipment in rear
- Wear PPE
- Slow hand & arm movements
- Do not block grilles
- Remove contaminated items after decontamination or place in sealed biohazard bags
- Do not store items in a biosafety cabinet
- Disinfect after completion & autoclave wastes



## BIOSAFETY CABINETS

- Use of flammable gases is not recommended in recirculating BSC's (Class II Type A1 and Type A2) because there could be a build up of gas if there is a leak or if the line is not completely turned off.
  - Several BSC have blown up.....
- Heat in a BSC & Safety
  - May cause turbulence within the air curtain which disrupts the pattern of HEPA-filtered air
  - May impact the integrity of the HEPA filter
  - May impact the seals of the cabinet



SURVEY OF GOVERNMENT AGENCIES & BSC  
MANUFACTURERS STRONGLY RECOMMEND  
AGAINST THE USE OF FLAMMABLE GASES AND  
FLAMES

- NIH/CDC
- WHO
- Public Health Agency of Canada
- NSF/ANSI (international standards)
- The Baker Company
- NuAire



## ALTERNATIVES TO PLUMBED-IN NATURAL GAS

- FirBoy
- Bacti-Cinerator IV
- Glass Bead Sterilizer
- Electric Bunsen Burner



## BIOSAFETY CABINETS....MORE INFORMATION

- Location of biosafety cabinet
- Annual certification
- Re-certify if moved, or see a drop in air pressure gauge
- Decontamination before moving/discarding
- Use of UV light:
  - Hazard to eyes/skin/cabinet
  - 1-2 year life
  - Poor Penetration-dust protects from deactivation
  - 10% bleach with a 70% ethanol wipe

